



## Is an Elbow Orthosis the Right Solution For You and Your Dog?

Elbow pathology is common. The most common reasons for considering an elbow orthosis include dynamic luxation (dislocation) and osteoarthritis.

Clinical signs of elbow injury depend on severity but may include lameness, swelling, and mal-alignment. Elbow luxation is usually managed with sedation and reduction +/- surgery. Chronic cases are sometimes managed with an orthosis. Elbow osteoarthritis is common and one of the most challenging joints to manage. There is no single effective therapy. Management must be by an integrative approach and include: weight management, medical management, supplements, rehabilitation, surgery, partial joint replacement, and sometimes an orthosis.

A custom orthosis can play a role in therapy and has the following purposes:

- 1) Restrict rotation at the elbow. The elbow is composed of 3 bones (humerus, radius and ulna). The radius and ulna compose the forearm and this allows some rotation around the humerus. This type of motion can induce pain in the arthritic joint.
- 2) Restrict frontal plane movement (varus and valgus). These are sideways movements. These movements are very subtle in the canine elbow and as such not as important as rotational movement can be in causing pain.
- 3) Limit joint motion (flexion and extension) to the most comfortable range. The device can be designed to allow movement within controlled parameters because end range flexion and end range extension can be painful.

Custom elbow orthoses cannot take away all pain associated with elbow pathology. Animals must bear weight through the affected joint. However, by limiting the above the device is helpful when your dog is active. In this way the device is a support device.

Because an orthosis is not the correct therapy for all patients, before choosing an orthosis the following points are important to keep in mind:

- 1. Device design is paramount to success.** Careful consideration is taken in prescribing a device and its specific components.
  - a. Articulation (hinging) is ideal whenever possible in order to provide as close to normal limb use as possible. Articulation can be controlled with design and the use of customized hinges and straps.
  - b. Not all dogs are candidates for elbow orthoses because of body shape. Dogs with a deep chest, short humerus (upper arm) or excessive weight are not candidates because the device cannot be fitted properly.



2. **The device MUST be put on and removed daily.** The orthosis stabilizes the elbow from the outside, while surgery does so from the inside. Therefore it is important to follow the exact wearing schedule provided by OrthoPets and A Loyal Companion. Wearing schedules vary with type of injury.
3. **Adjustments are expected and are a normal part of the custom orthosis process.** The device is custom-made for your pet. Every effort is made to accurately fit the device. If adjustments are necessary to optimize the fit and function of the device, A Loyal Companion can complete most adjustments in Tucson. Those devices requiring more extensive adjustment will be shipped to OrthoPets and require a 1-3 day turn around time, not including shipping time. Your pet is much more active at home than at the clinic. Increased activity can expose fit concerns requiring adjustment. Additional adjustments are most commonly required in the first few months based on follow up appointments with A Loyal Companion. Adjustments are a normal part of a custom orthosis. While adjustments are not always required, they should be anticipated. (see importance of follow-up #4). Please follow all instructions with regard to monitoring the leg and contact A Loyal Companion promptly if you have concerns.
3. **Follow-up is critical to success.** Proper use of the orthosis is necessary to meet therapeutic goals and to ensure its safe application over the lifetime of your dog or the duration of injury healing. In the first few months of fitting your doctor with the help of A Loyal Companion will see you and your dog for fit checks and coaching with regard to device use. Annual to twice annual appointments, depending on injury, age and activity of your dog, are needed. At these appointments your doctor will thoroughly assess your dog's orthopedic condition and evaluate the condition/fit of the device. Recommendations will be made for continued success in the device.
4. **Rehabilitation, the first key for success.** Most dogs adapt quickly to wearing an orthosis. Behavioral techniques can facilitate this. Also, your dog will need to learn basic skills while wearing the device. These include: transitions (sitting, lying down, and getting up), stairs, getting into vehicles safely, managing on different types of surfaces (ground, carpet, hardwood floor, etc.). Finally, orthopedic injury leads to compensatory abnormal movement and associated muscle strain and weakness. The best way to ensure the highest level of success with an orthosis is to follow a rehabilitation schedule. Each patient's condition and abilities are unique and as such an individualized rehabilitation program is needed. OrthoPets and A Loyal Companion strongly advise working with a certified canine rehabilitation professional (CCRT or CCRP). Please consult your family doctors or surgeon for referral to a veterinary rehabilitation professional in your area.
5. **A proactive approach to arthritis management is the second key to long-term success.** If the joint itself is injured rather than a ligament alone, osteoarthritis may develop. Just as rehabilitation is important, arthritis management is key as well. Steps



PH 520-448-3548  
FX 520-457-5828

4758 E. Grant Road  
Tucson, AZ 85712

info@aloyalcompanion.com  
www.aloyalcompanion.com

taken early and continued throughout your dog's lifetime will make a difference in terms of regaining and maintaining comfort and an active lifestyle well into the senior years. Consult with your doctors for a comprehensive pain management plan.



PH 520-448-3548  
FX 520-457-5828

4758 E. Grant Road  
Tucson, AZ 85712

info@aloyalcompanion.com  
www.aloyalcompanion.com



PH 520-448-3548  
FX 520-457-5828

4758 E. Grant Road  
Tucson, AZ 85712

[info@aloyalcompanion.com](mailto:info@aloyalcompanion.com)  
[www.aloyalcompanion.com](http://www.aloyalcompanion.com)